







### Guideline for 1:1 Drawing

1. Fill in all the measurements in (  mm  inch ) at the following pages that requested as below.
2. Mark sensitive areas with a "+" (plus sign) on the diagram.
3. The Trial Prosthesis fitting is most successful when the greatest number of measurements are recorded and may be able to make slight modifications to the order form.
4. All lengths should be measured from the palmar side.

 = Length

 = Circumference



This is an example, please draw the hand outline and write the measurements on the next pages.






Please draw 1:1 outline of the sound side (use separate sheet if necessary) and enter measurements as requested by the guideline on page 3.

# Palmar Side Face Down (Sound Side)

Example



\*This information does not replace cast. It should be provided as additional information.