



Semi-Custom Made (SCM) Order Form Model TOE Series



Company name, Email, Tel, Fax, Contact person / Account no., Order date, PO no.

Disclaimer - Please explain to the patient and ensure the patient understands that:

- 1. It is not possible to fabricate the appearance of the prosthesis exactly same as the sound side.
2. The size of the prosthesis may not match to the sound side, it depends on the condition of the residual limb.
3. The color of our silicone prosthesis may not match to the patient's skin tone and that, please therefore do not expect that the silicone prostheses match the patients' skin tone at all time.
4. The lives of the prosthesis depend on the environment, usage and maintenance.
5. The warranty policy of Regal silicone prosthesis, please refer to Regal catalogue 2015 page 154-155.
6. The leadtime is 14 - 21 working days for Trial Prosthesis, and 14 working days for Final Prosthesis. Delivery time may be affected by the missed or rescheduled appointments, delayed shipment or other causes beyond our control.
7. The Trial Prosthesis is strictly for the purpose of checking the color and the suitability, not for long term use or re-sell, and should be return to Regal upon requested.

Prosthetist Signature, Date


Ordering procedure:

- 1. Take applicable external measurements of the patient's sound and residual sides; draw 1:1 outline and enter data into the following pages.
2. Using the color sample from Regal, select the color that closely resembles the patient's sound and residual sides and enters into the following pages.
3. Take photos of the residual and sound sides.
4. Using the photos in the catalogue (2015), select the size that most closely resembles the patient's sound side and enters into the following pages.
5. Cast the sound and residual sides in the preferred natural position.
6. Mark casts with the patient name, APML alignment and sensitive area.
7. Securely wrap the cast before sending out. As casts can be broken easily during transportation.




### Guideline for 1:1 Drawing

1. Fill in all the measurements in (  mm  inch ) at the following pages that requested as below.
2. Mark sensitive areas with a "+" (plus sign) on the diagram.
3. The Trial Prosthesis fitting is most successful when the greatest number of measurements are recorded and may be able to make slight modifications to the order form.

 = Length

 = Circumference






Please draw 1:1 outline of the Residual and Sound sides (use separate sheet if necessary) and enter measurements as requested by the guideline on page 3.

# Plantar Side Face Down (Residual and Sound Sides)

Example



\*This information does not replace cast. It should be provided as additional information.